

TECHNICAL DATA SHEET

SANCARE Gynecological 450mm Natural PWD

STERILE LATEX GYNAECOLOGICAL GLOVES BEADED 450 mm

Powdered

Physical Dimensions

Size	Length (mm)	Palm Width (mm)	Thickness (mm) Minimum	
	Min		Palm	Finger
<i>Small</i>		80 ± 10		
<i>Medium</i>	450 mm	95 ± 10	0.15	0.15
<i>Large</i>		111 ± 10		
<i>Xtra Large</i>		114 ± 10		

Physical Properties

Characteristics	Before Ageing	After Ageing
Tensile Strength (Mpa) Min	14	14
Ultimate Elongation (%) Min	650	500
Force at Break (N) Min	6	6

Formulation

Curing Agent	<i>Sulphur</i>
Accelerators	<i>Zinc diethyl dithio carbamate, Zinc mercapto benzthiazole</i>
Activators	<i>Zinc Oxide</i>
Antioxidants	<i>Wingstay L (Condensation product of p-cresol & dicyclopentadiene</i>
Stabiliser	<i>Potassium Hydroxide</i>
Coagulant	<i>Calcium Nitrate</i>
Donning Powder	<i>Corn Starch (US FDA approved)</i>

Performance Requirements

Characteristics	Related Defects	Inspection Level	AQL
<i>Water tightness (Final Release)</i>	<i>Holes</i>	<i>G-1 Single Sampling Plan</i>	1.5

Primary Material	<i>Latex</i>
Colour	<i>Natural</i>
Shape	<i>Ambidextrous</i>
Cuff	<i>Beaded</i>
If Sterilised	<i>Ethylene Oxide</i>
Complies to	<i>EN 455 , ASTM 3578</i>



Packing

1 Pair in a Paper Wallet, then to a Paper Peel Down to Open Pouch

Quantity in 1 Inner Box	:	25 Pairs
Quantity in 1 Master Shipper Carton	:	100 Pairs (4 Inner Boxes)

Manufactured By : **SANREA Healthcare Products Pvt Limited**
KINFRA Integrated Industrial & Textile Park
Plot # P-56, Pearl Road, Kanjikode East PO,
Palakkad 678621, Kerala, INDIA
Tel : 91 491 2970145
Email : marketing@sanrea.com ; sanreamarketing@gmail.com
Web : www.sanrea.com

