

TECHNICAL DATA SHEET

SANCARE Gynecological 450mm Natural PF

STERILE LATEX GYNAECOLOGICAL GLOVES BEADED 450 mm

POWDER FREE CHLORINATED

Physical Dimensions

Size	Length (mm) Min	Palm Width (mm)	Thickness (mm) Minimum	
			Palm	Finger
Small		80 ± 10		
Medium	450 mm	95 ± 10	0.15	0.15
Large		111 ± 10		
Xtra Large		114 ± 10		

Physical Properties

Characteristics	Before Ageing	After Ageing
Tensile Strength (Mpa) Min	14	14
Ultimate Elongation (%) Min	650	500
Force at Break (N) Min	6	6

Powder Content

< 2mg /dm²

Protein Content

<50µg/dm²

Formulation

Curing Agent	Sulphur
Accelerators	Zinc diethyl dithio carbamate, Zinc mercapto benzthiazole
Activators	Zinc Oxide
Antioxidants	Wingstay L (Condensation product of p-cresol & dicyclopentadiene
Stabiliser	Potassium Hydroxide
Coagulant	Calcium Nitrate

Performance Requirements

Characteristics	Related Defects	Inspection Level	AQL
Water tightness (Final Release)	Holes	G-1 Single Sampling Plan	1.5

Primary Material	Latex
Colour	Pale Yellow
Shape	Ambidextrous
Cuff	Beaded
If Sterilised	Ethylene Oxide
Complies to	EN 455 , ASTM 3578



Packing

1 Pair in a Paper Wallet, then to a Paper Peel Down to Open Pouch

Quantity in 1 Inner Box	:	25 Pairs
Quantity in 1 Master Shipper Carton	:	100 Pairs (4 Inner Boxes)

Manufactured By : **SANREA Healthcare Products Pvt Limited**
KINFRA Integrated Industrial & Textile Park
Plot # P-56, Pearl Road, Kanjikode East PO,
Palakkad 678621, Kerala, INDIA
Tel : 91 491 2970145
Email :marketing@sanrea.com ; sanreamarketing@gmail.com
Web : www.sanrea.com