

TECHNICAL DATA SHEET

SANCARE Long Cuff Ambidextrous 450mm Natural PF

LATEX GLOVES, AMBIDEXTROUS LONG CUFF 450MM LENGTH

POWDER FREE CHLORINATED

Physical Dimensions

Size	Length (mm) Min	Palm Width (mm)	Thickness (mm) Minimum	
			Palm	Finger
<i>Small</i>		80 ± 10		
<i>Medium</i>	450 mm	95 ± 10	0.15	0.15
<i>Large</i>		111 ± 10		
<i>Xtra Large</i>		114 ± 10		

Physical Properties

Characteristics	Before Ageing	After Ageing
Tensile Strength (Mpa) Min	14	14
Ultimate Elongation (%) Min	650	500
Force at Break (N) Min	6	6

Powder Content

< 2mg /dm²

Protein Content

<50µg/dm²

Formulation

Curing Agent	<i>Sulphur</i>
Accelerators	<i>Zinc diethyl dithio carbamate, Zinc mercapto benzthiazole</i>
Activators	<i>Zinc Oxide</i>
Antioxidants	<i>Wingstay L (Condensation product of p-cresol & dicyclopentadiene</i>
Stabiliser	<i>Potassium Hydroxide</i>
Coagulant	<i>Calcium Nitrate</i>

Performance Requirements

Characteristics	Related Defects	Inspection Level	AQL
<i>Water tightness (Final Release)</i>	<i>Holes</i>	<i>G-1 Single Sampling Plan</i>	<i>1.5</i>

Primary Material	<i>Latex</i>
Colour	<i>Pale Yellow</i>
Shape	<i>Ambidextrous</i>
Cuff	<i>Beaded</i>
If Sterilised	<i>Ethylene Oxide</i>
Complies to	<i>EN 455 , ASTM 3578</i>



Manufactured By :

SANREA Healthcare Products Pvt Limited
KINFRA Integrated Industrial & Textile Park
Plot # P-56, Pearl Road, Kanjikode East PO,
Palakkad 678621, Kerala, INDIA
Tel : 91 491 2970145
Email :marketing@sanrea.com ; sanreamarketing@gmail.com
Web : www.sanrea.com