

# TECHNICAL DATA SHEET

## SANCARE Gynecological 400mm Natural PWD

STERILE LATEX GYNECOLOGICAL GLOVES BEADED 400 mm

PRE POWDERED

### Physical Dimensions

Size	Length ( mm ) Min	Palm Width	Thickness ( mm ) Minimum		
			Cuff	Palm	Finger
7.0		89 ± 6			
7.5	400	95 ± 6	0.13	0.17	0.19
8.0		102 ± 6			
8.5		108 ± 6			

### Physical Properties

Characteristics	Before Ageing	After Ageing
Tensile Strength ( Mpa ) Min	24	18
Ultimate Elongation (%) Min	750	560
Stress at 500% elongation (Mpa) Max	5.5	NA
Force at Break (N) Min	9	9

### Powder Content

<15 mg /dm<sup>2</sup>

### Protein Content

<200 µg/dm<sup>2</sup>

### Formulation

Curing Agent	Sulphur
Accelerators	Zinc diethyl dithio carbamate, Zinc mercapto benzthiazole
Activators	Zinc Oxide
Antioxidants	Wingstay L ( Condensation product of p-cresol & dicyclopentadiene
Stabiliser	Potassium Hydroxide
Mould Cleaner	Nitric Acid
Coagulant	Calcium Nitrate
Donning Powder	Corn Starch ( US FDA Approved )

### Performance Requirements

Characteristics	Related Defects	Inspection Level	AQL
Water tightness ( Final Release)	Holes	G-1 Single Sampling Plan	1.5

Primary Material	Latex
Colour	Natural
Shape	Anatomic
Cuff	Beaded
If Sterilised	Ethylene Oxide
Complies to	EN 455 , ASTM 3577



### Packing

1 Pair in a Paper Wallet, then to a Paper Peel Down to Open Pouch

Quantity in 1 Inner Box	:	50 Pairs
Quantity in 1 Master Shipper Carton	:	200 Pairs ( 4 Inner Boxes )

Manufactured By : **SANREA Healthcare Products Pvt Limited**  
KINFRA Integrated Industrial & Textile Park  
Plot # P-56, Pearl Road, Kanjikode East PO,  
Palakkad 678621, Kerala, INDIA  
Tel : 91 491 2970145  
Email : marketing@sanrea.com ; sanreamarketing@gmail.com  
Web : www.sanrea.com