

# TECHNICAL DATA SHEET

## SANCARE Surgical Natural Powder FREE

STERILE LATEX SURGICAL GLOVES BEADED

POWDER FREE CHLORINATED

### Physical Dimensions

Size	Length ( mm )	Palm Width ( mm )	Thickness ( mm ) Minimum		
	Min		Cuff	Palm	Finger
6.0		76 ± 6			
6.5		83 ± 6			
7.0		89 ± 6			
7.5	280	95 ± 6	0.10	0.10	0.10
8.0		102 ± 6			
8.5		108 ± 6			
9.0		114 ± 6			

### Physical Properties

Characteristics	Before Ageing	After Ageing ( 70° C for 166 Hrs )
	Tensile Strength ( Mpa ) Min	24
Ultimate Elongation (%) Min	750	560
Stress at 500% elongation (Mpa) Max	3.5	NA
Force at Break (N) Min	9	9

**Powder Content** <2 mg /dm<sup>2</sup>

**Protein Content** <50 µg/dm<sup>2</sup>

### Formulation

Curing Agent	Sulphur
Accelerators	Zinc diethyl dithio carbamate, Zinc mercapto benzthiazole
Activators	Zinc Oxide
Antioxidants	Wingstay L ( Condensation product of p-cresol & dicyclopentadiene
Stabiliser	Pottassium Hydroxide
Coagulant	Calcium Nitrate Corn Starch ( US FDA Approved )

### Performance Requirements

Characteristics	Related Defects	Inspection Level	AQL
Water tightness ( Final Release )	Holes	G-1 Single Sampling Plan	1.5

<b>Primary Material</b>	Latex
<b>Color</b>	Pale Yellow
<b>Shape</b>	Anatomic
<b>Cuff</b>	Beaded
<b>If Sterilized</b>	Ethylene Oxide / GAMMA
<b>Complies to</b>	EN 455 , ASTM 3577



### Packing

1 Pair in a Paper Wallet, then to a Paper Pouch

Quantity in 1 Inner Box	:	50 Pairs	
Quantity in 1 Master Shipper Carton	:	400 Pairs	( 8 Inner Boxes )

Manufactured By :

**SANREA Healthcare Products Pvt Limited**  
KINFRA Integrated Industrial & Textile Park  
Plot # P-56, Pearl Road, Kanjikode East PO,  
Palakkad 678621, Kerala, INDIA  
Tel : 91 491 2970145  
Email : marketing@sanrea.com ; sanreamarketing@gmail.com  
Web : www.sanrea.com